





Using your Handihaler®




1 **Open** the Handihaler® by pulling up on the dust cap. To expose the chamber where you will place your medication capsule, pull up on the mouthpiece.




2 **Place** the capsule in the center of the chamber.




3 **Close** the mouthpiece until you hear a "click".



4 **Holding** upright, push the green button on the side of the Handihaler® just once. This will break the capsule of medication.



5 **Breathe out** completely, away from the mouthpiece.



6 **Seal** your lips around the mouthpiece. **Breathe in fast and deep.** You should hear the capsule vibrate.

- ! 24/7/365 Emergency Support
- ☎ Toll Free: 1-855-672-6262
- ➔ CareicaHealth.com/oxygen



7




hold your
breath

Remove the HandiHaler® from your mouth. Hold your breath for 5-10 seconds.

You can repeat steps 5 to 7 to make sure you got all the medication.


8



Pull up on the mouthpiece and throw away the used capsule.


DO NOT touch the capsule. If you do, wash your hands; **DO NOT** touch your eyes.

9



Close the mouthpiece and dust cap.

10



Rinse your mouth with water; spit the water out. (If your medication is a corticosteroid e.g. flovent).

Tip

Clean by removing the capsule and rinsing with warm tap water. Allow to air dry.



24/7/365 Emergency Support



Toll Free: 1-855-672-6262



CareicaHealth.com/oxygen

